HQ-AT20

**MANUAL (p. 2)**  
BREATH ALCOHOL TESTER

**ANLEITUNG (S. 8)**  
ATEMALKOHOLTESTER

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
**MANUAL DE UTILIZARE (p. 62)**  
TESTER NIVEL ALCOOLEMIE

## ENGLISH


- |                         |                 |
|-------------------------|-----------------|
| (1) Power Button        | (2) LCD display |
| (3) Torch Light         | (4) Breath Pipe |
| (5) Torch Switch        | (6) Key chain   |
| (7) Battery Compartment |                 |

### Features

- 3-step alcohol test with audible alert
- LCD display
- Quick response and resume
- Auto power off
- Low battery indication
- Torch function
- Key chain
- Use 2 x AAA alkaline batteries

### Battery installation

- Slide open the battery compartment lid from the back of the unit.


- Insert two AAA (or UM-4) size alkaline batteries with correct polarity as indicated.
- Replace the battery compartment lid.

## How to use the alcohol tester

### Warm up

1. Press and hold the power button for 1 second, the LCD display will turn on with one beep sound.

2. The @ready@ symbol will flash for 15-20 seconds that indicate

the tester is in warm-up mode.


## DISPLAY ALCOHOL BREATH TESTER

### Note:

- \* If warm-up is defeated, the tester will be turned off automatically, you can repress the power button and try again.
- \* It always has a large number of gas & impurity absorbed on the surface of the alcohol sensor, when taking the tester out of the package for the first time, or if it hasn't been in use for a long while, you need to clean and warm the sensor up.

### Test alcohol concentration

1. When warm-up is succeed, the Test... symbol will be displayed that indicate the tester is in testing mode.

2. Blow into the breath pipe for 3-5 seconds continuous.
3. The test result is viewable on the LCD display.

4. The LCD keep the primary format (Safety symbol) when the alcohol concentration is the level:  
 Under 0.02% BAC (0.1mg/l BRAC)  
 Under 0.05% BAC (0.25mg/l BRAC)

5. The Caution symbol is displayed when the alcohol Concentration is the level:  
 0.02% to 0.05% BAC (0.1mg/l to 0.25mg/l BRAC)  
 0.05% to 0.08% BAC (0.25mg/l to 0.4mg/l BRAC)


### **Precaution**

- If you have consumed alcohol, it is recommended to wait with the test for 20 minutes. It takes the body approximately 20 minutes to absorb the alcohol into the blood from the digestive system, and residues of alcohol from the mouth takes this long to dissipate.
  - For an accurate result, do not repeat next testing within 3 minutes after last time.
  - Do not smoke and eat for 20 minutes prior testing.
  - Do not use chemicals to clean the tester directly.
  - Do not blow liquid into the breath pipe.
  - Do not place the tester in the closed environment with the contaminated air.
  - The result of this alcohol gives BAC for reference only, manufacturer does not to take any legal responsibility.
6. The danger symbol is displayed with beep alert when the alcohol concentration has the level:  
 Over 0.05% BAC (0.25mg/l BRAC)  
 Over 0.08% BAC (0.4mg/l BRAC)

### Low battery indication

If (low battery) indicator show on the display, please try to replace the old batteries (use the alkaline battery only).

### Auto power off

1. Power will automatically turn off after the test lasts for about 25 seconds.
2. If warm-up is defeated, the tester will automatically turn off also.

### **How to use the torch**

Set the torch switch to <on> position for turn on the torch.

- \* Avoid inaccurate result, do not test the alcohol concentration when the torch is lit.

### **Specifications**

- \* Sensor: Advanced semiconductor oxide alcohol sensor
- \* Range: 0.00%-0.1% BAC (0.00-0.50mg/l BRAC)
- \* Accuracy: 0.01% BAC (0.05mg/l BRAC)
- \* Power input: 2 X 1.5V @AAA@ alkaline battery
- \* Operating time: warm up: 15-20s, response: <5s
- \* Auto power off: 1 minute approx.
- \* Dimension: 100X45X20mm

**DECLARATION OF CONFIRMITY**

We,

Nedis B.V.

De Tweeling 28

5215MC 's-Hertogenbosch

The Netherlands

Tel.: 0031 73 599 1055

Email: info@nedis.com

Declare that product:

Name: HQ

Model: HQ-AT20

Description: Breath alcohol tester with LCD display

Is in conformity with the following standards

EMC: EN 55022: 1998 + A1: 2000 + A2: 2003

EN 61000-6-3: 2001

EN 61000-6-1: 2001

Following the provisions of the 89/336/EEC directive.

's-Hertogenbosch, 2008, February 26

Mrs. J. Gilad

Purchase Director

**Safety precautions:**

To reduce risk of electric shock, this product should **ONLY** be opened by an authorized technician when service is required. Disconnect the product from mains and other equipment if a problem should occur. Do not expose the product to water or moisture.

**Maintenance:**

Clean only with a dry cloth. Do not use cleaning solvents or abrasives.

**Warranty:**

No guarantee or liability can be accepted for any changes and modifications of the product or damage caused due to incorrect use of this product.

**General:**

Designs and specifications are subject to change without notice. All logos brands and product names are trademarks or registered trademarks of their respective holders and are hereby recognized as such.


**Attention:**

This product is marked with this symbol. It means that used electrical and electronic products should not be mixed with general household waste.

There is a separate collections system for these products.

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Copyright ©
Copyright © (Autorská práva)