**KONIG**  
HEALTH CARE

**HC-PS500**

**Personal scale**

**MANUAL**

**ANLEITUNG**

**MODE D'EMPLOI**

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**MANUALE**

**MANUAL DE USO**

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**KÄYTTÖOHJE**

**BRUKSANVISNING**

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**MANUAL DE UTILIZARE**

English

Deutsch

Français

Nederlands

Italiano

Español

Magyar

Suomi

Svenska

Česky

Română

**Advice for use and care:**

Ensure the scale is on a flat, even and firm surface. Avoid carpet or soft surfaces. Always weigh yourself with the same scale placed on the same surface. Uneven floors may affect the reading.

The surface of the scale will be slippery if wet. Keep it dry!

Weigh yourself without clothing and footwear, before meals and always at the same time of day. Keep the scale in a dry place to protect the electronic components. Clean after use with a lightly damp cloth, do not use solvents or immerse the product in water. Avoid overloading the scale, otherwise, permanent damage may occur.

A user who has any metal material inserted into his/her body may get differences in measurement of body fat reading. Pacemaker users are advised against using this device.

The scale is designed for household use only, not to be used for commercial purposes.

**Weighing operation:**

1. Remove the insulation tab from the battery contacts before use.
2. Select the weighing mode using the switch (Kg/St/Lb) located on the back of the scale.
3. Step on the scale platform to switch on the product, stand still while the scale computes your weight.
4. Your weight will be displayed and fixed.
5. When you step off the scale it will automatically switch off after 6 seconds.
6. If no weighing is carried out while 0.0 is displayed the scale will automatically switch off after 6 seconds.

Remarks: This initialization process as below must be repeated if the scale is moved. At all other times step straight on the scale.

- Step 1: Press the platform centre and remove your foot.
- Step 2: "0.0" will be displayed
- Step 3: The scale will switch off and now it is ready for use

**Body fat measurement:**

In order to ensure accurate results users should enter the required personal parameters: age, sex and height.

These scale has the facility to store the personal parameters of 10 people.

The technology used in these scale, bio-electrical impedance analysis (BIA) is more technologically advanced than the older Body Mass Index (BMI) method where fat % is calculated using only height and weight. BIA technology uses 5 parameters to get a more accurate resolution for the analysis (weight and biological resistance which are automatically measured by the scale together with the three personal parameters entered manually).

Mode selection and parameter setting:

Press and hold **O** to switch on/ switch off the scale.

Press **▲** or **▼** to switch between body parameter 0-9. Press **■** to confirm.

Press **▲** or **▼** to select male or female. Press **■** to confirm.

Press **▲** or **▼** to indicate height. Press **■** to confirm.

Press **▲** or **▼** to indicate age. Press **■** to confirm.

The values are flashed after 0.0 appears. First the weight is recorded, then  is displayed, indicating that the scales are analyzing, please do not alight. Once finished the weight and body fat & water will be displayed alternately for 4 times before the scales automatically turn off. Compare the body fat reading against the (body fat & water ratio chart) to confirm the body status.

Body fat ratio chart

AGE	FEMAL (%)				MALE (%)			
	Low	Healthy	High	Obese	Low	Healthy	High	Obese
10-39	<21	21-33	33-39	>39	<8	8-19	19-25	>25
40-59	<23	23-35	35-40	>40	<11	11-22	22-28	>28
60-79	<24	24-36	36-42	>42	<13	13-25	25-30	>30

Due to the natural differences between people this chart should be used for reference purposes only.

Water ratio chart & muscle ratio chart (%)

AGE	FEMAL (%)				MALE (%)			
	Low	Healthy	High	Body Muscle	Low	Healthy	High	Body Muscle
10-15	<57	57-67	>67	>34	<58	58-72	>72	>40
16-30	<47	47-57	>57		<53	53-67	>67	
31-60	<42	42-52	>52		<47	47-61	>62	
61-80	<37	37-47	>47		<42	42-56	>56	

Bone weight chart (KG)

	femal			male		
weight	<45Kg	45~60Kg	>60Kg	<60Kg	60~75Kg	>75Kg
Ref. bone weight	~5,8Kg	~6,2Kg	~6,5Kg	~6,5Kg	~6,9Kg	~7,2g

**Specifications:**

Equipped with 4 precision load cells

Digital LCD display

Capacity 180Kg/396lb

Graduation: 0.1 Kg / 0.2lb / 1/4lb

Body fat division: 0.1%

Water division: 0.1%

Age range for body fat: 10 - 80 years

Height range: 100 - 220cm / 3ft 3.5in - 7ft 2.5in

10 person memory

Choice of weighing modes: Kilogram (kg), stone (st), Pound (lb)

Operates with 3 x 1.5V AAA battery

**Declaration of conformity**

We,  
Nedis B.V.  
De Tweeling 28  
5215MC 's-Hertogenbosch  
The Netherlands  
Tel.: 0031 73 599 1055  
Email: info@nedis.com

Declare that product:

Name: KONIG  
Model: HC-PS500  
Description: Electronic personal scale

Is in conformity with the following standards

EMC: EN 61000-6-3 : 2001 + A11 : 2004  
EN 61000-6-1 : 2001

Following the provisions of the 2004/108/EC directive.

's-Hertogenbosch, 2007, November 16

Mrs. J. Gilad  
Purchase Director

**Safety precautions:**

This product should **ONLY** be opened by an authorized technician when service is required. Do not expose the product to water or moisture.

**Maintenance:**

Clean only with a dry cloth. Do not use cleaning solvents or abrasives.

**Warranty:**

No guarantee or liability can be accepted for any changes and modifications of the product or damage caused due to incorrect use of this product.

**General:**

Designs and specifications are subject to change without notice.

All logos brands and product names are trademarks or registered trademarks of their respective holders and are hereby recognized as such.

**Attention:**

This product is marked with this symbol. It means that used electrical and electronic products should not be mixed with general household waste.

There is a separate collections system for these products.

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Copyright ©